

## Infusion stand

Both our standard infusion stand and our dosage pump infusion stand are designed in keeping with the qualitative and functional standards demanded by Swedish hospital professionals.

### Frame

All-welded chromium-plated steel tubes with four hooks for intravenous containers.

### Locking device

The stand is locked at the desired height with an expansion-locking device.

The expander clamps firmly around the whole inner tube, which ensures safe locking and also extends the useful life of the stand.

### Base

The heavy; detachable, five-leg wheel base provides stability.

### Wheels

Easy-roll castors glide smoothly and safely also on carpeted areas.

### Steering handle

Vertically adjustable and detachable.

### Infusion stands for dosage pump

Has a stable construction with the dosage pump mounted at the centre of the stand.

### Dimensions and specifications

Heights of standard infusion stand: 125-225 cm.

Heights of infusion stand for dosage pump: 129-229 cm.

Diameter of wheel base: 59 cm.

Wheels: 100-mm castors with grey rubber surface. Antistatic castors as accessories.



Infusion Stand	Item.No.
Standard infusion stand with steering handle	080-3110
Standard infusion stand without steering handle	080-3050
Infusion stand for dosage pump with steering handle	080-3200
Infusion stand for dosage pump without steering handle	080-3210
<b>Accessories</b>	
Antistatic castors	080-3001
Steering handle	080-3002

THE INFUSION STANDS are delivered in simple assemble packs to reduce transportation overheads/delivery costs.

*Subject to alternations without notice.*